



tangible™
HYDRA-PEG

Erfahrungsbericht 12 Monate

Michael Wyss


M.Sc. Optometrist FAAO

eyeness ag, Bern / Switzerland

mwyss@eyeness.ch




Tangible™ Hydra-PEG



Tangible™ Hydra-PEG

Tangible Hydra-PEG Coating




OPTIMUM | tangible™
HYDRA-PEG

- Polymer encapsulates lens
- Improves wettability
- Increases surface water retention
- Increases lubricity
- Reduces deposits
- Improves comfort

Resultat 14 Tage


HydraPEG Contact Lens Surfacing on Oasys (JNJ)

P. Caroline, et al., Pacific University




HydraPEG Contact Lens Surfacing on Optimum Extra

M. Wyss et al., eyeness ag N = 34 eyes




Resultat 14 Tage

34 Augen – nach Patienten Typen aufgeteilt




Resultat 14 Tage

Case Report EJ, Keratokonus ohne Hydra-PEG



Resultat 14 Tage

Case Report EJ, Keratokonus mit Hydra-PEG






Tangible™ Hydra-PEG Coating Resists Lens Deposits

All GP Scleral Lenses Dipped in Lipid/Dye Mixture

Resultate 12 Monate


18 Augen




Haltbarkeit

Case Report EJ, Keratokonus


untreated



2 week Hydra-PEG



1 year Hydra-PEG



Pflegehinweise (Falco)



MULTIPURPOSE SOLUTIONS:

Boston Simplus®
Regard® K
Menicare Plus
Solocare Aqua
Perfect® cleaner / storage solution for
GP lenses

PEROXYD SYSTEME:

TAB IN ONE Peroxide System
Bilopur

ENZYMATISCHE REINIGER:

Concare Enzyme Tablets



ALKOHOLHALTIGE REINIGER:

Oté Clean
Contopharma i-clean! Alcohol Cleaner
Miraflow Cleaner und vergleichbare
Produkte

ABRASIVE REINIGER:

Boston Advance

DIVERSE:

Progent
Oté Twins Active
Boston One Step Liquid Enzymatic
Boston Lab Cleaner

TANGIBLE
HYDRA-PEG
BESCHICHTUNG

falco-linsen.com

falco

eyeness
wellness für ihre augen



www.eyeness.ch/news/downloads